


Digital Rotator CLR-01



Digital Rotator

Rotor is provides with
2 different types of
tube holders



2nd option-with disc tube holder

Features :

- ▶ Rotor is provided with 2 different types of tube holders.
- ▶ In 1st option rotor is supplied with rotisserie holder for tubes.
- ▶ In 2nd option rotor is supplied with disc holder for tubes & the disc angle can be set from 0-90degree.
- ▶ Different options of holding accessories.
- ▶ Speed & time adjustment.
- ▶ LCD with backlite displays time and speed.
- ▶ DC motor for reliable operation.
- ▶ Designed for continuous or time based operation.

Specifications :

- ▶ Tilt angle : 90degree
- ▶ Speed Range : 10-70rpm
- ▶ Timer : 1~1200min.
- ▶ Dimension (mm) : 195 (W) x 190 (L) x 260 (H)
- ▶ Weight : 8kg
- ▶ Power supply : 100~250 V AC, 50/60 Hz, 40 Watt

1st option-with rotisserie tube holder

Dimension (mm) : 205 (W) x 515 (L) x 280 (H)

2nd option-with disc tube holder

Dimension (mm) : 195 (W) x 270 (L) x 360 (H)

All Specifications and designs are subject to change for improvement.



CONTECH[®]
Instruments Ltd.
AN ISO 9001,14001 & OHSAS 18001 COMPANY
Where innovation never stops.....

301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.
Tel. :+91 22 6139 3000 (100 lines)
Fax :+91 22 2761 8421
E-mail : info@contechindia.in / sales@contechindia.in
Website: www.contechindia.com



State Of The Art Manufacturing Plant