

Digital Tube Roller CTR-01






Features :

- ▶ Provides gentle rocking & rolling motion.
- ▶ Rocking & rolling motion for complete mixing.
- ▶ Six no. easy cleaning rollers.
- ▶ Design for continuous or time based operation.
- ▶ DC motor for reliable operation.

Specifications :

- ▶ Amplitude : 24mm
- ▶ Maximum Load Capacity : 4kg
- ▶ Speed Range : 10-70rpm
- ▶ Time setting range : 1~1200min.

- ▶ Dimension (mm) : 250 (W) x 460 (L) x 110 (H)
- ▶ Weight : 5kg
- ▶ Power supply : 100~250 V AC, 50/60 Hz, 30 Watt



All Specifications and designs are subject to change for improvement.