


ANTI VIBRATION TABLE FOR ANALYTICAL BALANCE

1. Anti vibration table is ideal for R & D, testing & QC laboratories located in a setup subjected to high vibrations due to operation of machines.
2. Anti vibration table is designed for placing analytical & lab balances used for high precision weighing.
3. Anti vibration table top is provided with granite surface plate with high accuracy finish.
4. Damping of vibration is provided at two locations, one below the granite top surface & second at the bottom leveling feet which can also help in proper leveling of the table.



CONTECH[®]
Instruments Ltd.

AN ISO 9001,14001 & 45001 COMPANY

Where innovation never stops.....

301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.

Tel. :+91 22 6139 3000 (100 lines)

E-mail : sales@contechindia.in

Website: www.contechindia.com

FACTORY / NABL LAB



Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone),
Mhape, Navi Mumbai - 400 710.