



# Multi Station Magnetic Stirrer with Heating CMS-10H



## Features :

- ▶ High performance multi station magnetic stirrer with heating.
- ▶ 10 station position with uniform temperature distribution.
- ▶ Microprocessor speed controller guarantees speed accuracy.
- ▶ Backlite LCD display for temperature & speed indication.
- ▶ Brushless DC motor for maintenance free & noiseless operations.
- ▶ Inbuilt temperature safety cutoff.
- ▶ Stainless steel anti corrosion top.
- ▶ Work Plate Dimension : 200 x 470 (mm)
- ▶ Maximum stirring Quantity : 400ml in each position.
- ▶ Maximum Stirring Speed : 1100rpm
- ▶ Temperature range of work plate : Ambient to 120°C  
(Actual liquid temperature is lesser than the work plate temperature)

## Specifications :

- ▶ Dimension (mm) : 205(W) X 610(L) X 65 (H),
- ▶ Weight : 4Kg approx.
- ▶ Heater Power : 480 W
- ▶ Power supply : 220~240 V AC, 50/60 Hz,  
500 Watt

All Specifications and designs are subject to change for improvement.